

Standard Operation Procedure for BD LSRFortessa X-20

1. Check the sheath cube and the waste tank, and the FACS Flow Cart is on.

- If the sheath cube is low, remove the probe, change a new cube, replace the probe, and press “Restart” on the FACS Flow Cart.
- If the waste tank is full, remove the probe and empty the waste. Add 300 ml of bleach, replace the probe and press “Restart” on the FACS Flow Cart.



- ## 2. Turn on the computer. Log in to Windows as the BD operator. The password is “Welcome#1”. Double-click the DIVA icon to open the instrument acquisition software. Log in to the software.




- ## 3. Turn on the instrument system power by pressing the green button on the right side of the instrument.



4. Rinse the System (First user of the day)


- a) Install a tube containing 3 mL of 100% Contrad70 solution on the SIP.
- b) Press RUN and HI on the cytometer fluidics control panel.
- c) Allow the solution to run on HI for 10 minutes.
- d) After 10 minutes, change a tube containing 3 mL of BD Rinse buffer on the SIP.
- e) Allow the solution to RUN on HI for 10 minutes.
- f) After 10 minutes, change a tube containing 3 mL of DI water on the SIP.
- g) Allow the solution to RUN on HI for 10 minutes.
- h) Press STANDBY button on the fluidics control panel.
- i) Place a tube containing no more than 1 mL of DI water on the SIP.
- j) Allow 30 minutes for the lasers to warm up before running samples or CS&T.



5. Running CST beads (First user of the day) (20 points)

- a) Prepare CS&T research beads:
 - 350uL filtered PBS + 1 drop of CST bead
 - Vortex well

b) Select **Cytometer > CST**



- c) **Verify that the lot ID under Setup Beads** matches the information on the BD FACSDiva CS&T research beads vial.
- d) Install the pre-prepared CST bead tube on the cytometer
- e) Select the **Load Tube Manually** checkbox
- f) Remove the DI water tube from the SIP
- g) Set the flow rate to LO and press RUN on the cytometer fluidics control panel
- h) Click **Run** in the CST workspace, and click **OK** to confirm that the tube has been installed.
- i) When the performance check is complete, remove the CST tube from the cytometer.
- j) Install a tube with 1 mL DI water onto the cytometer and press STANDBY
- k) In the **Cytometer Setup and Tracking performance** check completion dialog, click **View Report**.
- l) Verify that the cytometer performance passed
- m) Select **File > Exit** to close the Cytometer Setup and Tracking workspace and reconnect to the BD FACSDiva interface.

Important Notes:

- **QC is run by Core staff.**
- **You only need to run the QC only if Core Staff have not run it for the day.**
- **If you run the QC and if CST fails please contact Core staff via email.**

6. Shut Down

- a) Install a tube containing 3 mL of 100% Contrad solution onto the cytometer, RUN HIGH for 5 minutes.
- b) Install a tube containing 3 mL of BD FACSClean solution onto the cytometer with the support arm to the side (vacuum on).
- c) Allow the solution to RUN on HI for 1 minute
- d) Move the support arm under the tube (vacuum off) and continue to RUN on HI for 5 minutes.
- e) Repeat b)-d) steps with full-strength BD FACSRinse solution
- f) Repeat step b)-d) steps with DI water.
- g) Install a tube with no more than 1 mL of DI water onto the cytometer
- h) Set the fluidics mode to STANDBY.
- i) Exit BD FACSDiva software.
- j) Turn off the cytometer.
- k) Log out the computer.

Important Notes:

- **DO NOT leave more than 1 mL of water on the SIP.** When the instrument is turned off or left in standby mode, a small amount of fluid will drip back into the sample tube. If there is too much fluid in the tube, it could overflow and affect the cytometer performance.
- **DO NOT turn off the power on the BD FACSFlow supply system (FFSS).**
- **Check the waste tank after use. Empty the waste tank if necessary.**